


# ASSAY PERFORMANCE SUMMARY

## Duloxetine



Method	LC/MS/MS
Matrix	human plasma
Anti-coagulant	K <sub>3</sub> EDTA
LLQ	0.200 ng/mL
ULQ	100 ng/mL
Sample Volume	200 µL
Internal Standard	Stable label
Accuracy: %Bias	
LLQ	-2.5
LQC	-6.0
MQC	-5.0
HQC	-5.4
DilQC (DF=10)	2.0
Precision: %CV	
LLQ	6.7
LQC	7.9
MQC	9.1
HQC	6.6
DilQC (DF=10)	3.0
Stability: Minimum	
Ambient	7 hours
-20°C	In Progress
No. Freeze/thaws	4