


ASSAY PERFORMANCE SUMMARY

Fluconazole



Method	LC/MS/MS
Matrix	human plasma
Anti-coagulant	K ₃ EDTA
LLQ	25.0 ng/mL
ULQ	10,000 ng/mL
Sample Volume	50.0 µL
Internal Standard	Stable label
Accuracy: %Bias	
LLQ	4.0
LQC	-0.7
MQC	7.3
HQC	1.6
DilQC (DF=10)	-7.4
Precision: %CV	
LLQ	6.7
LQC	4.6
MQC	3.5
HQC	3.1
DilQC (DF=10)	10.6
Stability: Minimum	
Ambient	48 hours
-20°C	133 days
No. Freeze/thaws	4