


ASSAY PERFORMANCE SUMMARY

Lopinavir



Method	LC/MS/MS
Matrix	human plasma
Anti-coagulant	K ₃ EDTA
LLQ	5.00 ng/mL
ULQ	5000 ng/mL
Sample Volume	100 µL
Internal Standard	Stable label
Accuracy: %Bias	
LLQ	-1.8
LQC	-1.3
MQC	1.0
HQC	5.0
DilQC (DF=10)	5.0
Precision: %CV	
LLQ	7.4
LQC	10.4
MQC	3.5
HQC	3.5
DilQC (DF=10)	1.7
Stability: Minimum	
Ambient	7 hours
-20°C	24 days
No. Freeze/thaws	4